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16 REFERENCIA

Cao X. COVID-19: immunopathology and its implications for therapy. Nat Rev Immunol. 2020;20:269–70.

STORAGE: Refrigeration at 2-8°C
EXPIRATION DATE: The date is indicated on the packaging box

TYPE OF SAMPLE: Serum
DIAGNOSTIC KIT FOR USE IN VITRO

CONTENTS	
1 INTENDED USE	8 SAFETY CONSIDERATIONS
2 USERS	9 SAMPLES
3 INTRODUCTION	10 ASSAY PROCEDURE
4 PRINCIPLE	10.1 PREPARATION OF REGENTS
5 KIT COMPONENTS	10.2 PREPARATION OF SAMPLE
6 REQUIRED MATERIALS NOT INCLUDED IN THE KIT	10.3 ASSAY PERFORMANCE
7 REAGENT PREPARATION AND STORAGE	11 RECOMMENDATIONS FOR THE ASSAY
7.1 WASHING SOLUTION	12 QUALITY CONTROL
7.2 DILUENT AGENT	13 RESULTS INTERPRETATION
7.3 DILUENT SOLUTION OF ANTIBODIES	14 TEST LIMITATIONS
7.4 SECONDARY ANTIBODY	15 DEVELOPMENT FEATURES
7.5 DEVELOPMENT SUSTRATES	15.1 SENSITIVITY
7.6 STOP SOLUTION	15.2 ACCURACY
7.7 POSITIVE CONTROL	15.3 CROSS REACTIVITY
7.8 NEGATIVE CONTROL	16 REFERENCE

Interpretation of the positive test for anti-SARS-CoV-2 antibodies
The presence of IgG-type antibodies suggests that a subject has been exposed to the virus and has developed an immune response, typically this occurs at least two weeks after exposure and clinical symptoms expression of the disease. It does not categorically determine that such subject is no longer at risk of disease, but it suggests that is at lower risk than someone lacking of antibodies.

1 INTENDED USE
The kit allows detection of IgG-type antibodies specific to receptor-binding domain (RBD) of the spike protein (S) of the Severe Acute Respiratory Syndrome-coronavirus type-2 (SARS-CoV-2) in serum samples of subjects suspected of had been infected by the virus. This assay only provides qualitative results. The result should be correlated by health personnel with the patient's clinical data and results of complementary studies.

SARS-CoV-2 (formerly known as 2019-nCoV) has been reported to infect human respiratory epithelial cells through interaction with human angiotensin-converting enzyme-2 (ACE2) receptor. Protein S is a type-I transmembrane protein, consisting of two subunits, S1 and S2. The S1 subunit contains the receptor-binding domain (RBD), which is responsible for the recognition of ACE2 on cell surface, while subunit S2 contains basic elements necessary for membrane fusion. Protein S plays a key role in inducing immune response of T cells and generating neutralizing antibodies, which provides a protective immunity. The methodology developed for this kit aims to detect the presence of IgG- type antibodies anti-RBD of SARS-CoV-2 coronavirus S1 protein by means of an indirect ELISA assay. For this purpose, an immunoassay plate is sensitized with the RBD of SARS-CoV-2 coronavirus S1 protein (antigen), and then serum samples are added to determine presence of antibodies recognizing specifically RBD. The interaction of protein with IgG-type antibodies in the tested serum is evidenced by addition of an anti-human-IgG antibody coupled to peroxidase and subsequent enzymatic reaction with chromogenic substrate detected at 450 nm.

2 USERS
The kit is for use by clinical laboratory and/or healthcare professionals.

3 INTRODUCTION
The COVID-19 disease (by its abbreviation Coronavirus Disease 2019) is an acute respiratory disease caused by SARS-CoV-2, and such was declared a pandemic by the World Health Organization (WHO) in early 2020.

Glycoprotein S of coronaviruses is essential for binding the virus to the host cell during the infection process.

4 PRINCIPLE
The present ELISA kit is designed for qualitative detection of serum human antibodies of anti-SARS-CoV-1 IgG type. The assay uses a plate technique based on an enzymatic-immunoassay. The tested serums are diluted 1:100 and placed on a 96-well immunoassay plate previously sensitized with the RBD of SARS-CoV-2 virus (antigen) S-protein, the non-binding proteins are removed by subsequent washings. Formation of anti-RBD antigen-antibody complex is showed by adding anti-human IgG secondary antibody coupled to horseradish peroxidase (HRP), excess of secondary antibody is removed by washings. Adding 3,3',5,5'-tetramethylbenzidine (TMB) chromogenic substrate with hydrogen peroxide shows formation of anti-human IgG antibody anti-SARS-CoV-2-RBD-IgG antibody complex, through color development which can be quantified using a microplate reader at 450 nm.

The amount of anti-RBD antibodies in the sample will be proportional to the optical density (OD) shown by the chromogenic substrate color.

5 KIT COMPONENTS (One plate version)

ITEM No.	COMPONENT	QUANTITY	UNITS	CONTENT	STORAGE TEMPERATURE
B-2020	Sensitized immunoassay plate	1	1 plate	Nunc-MaxiSorp 96-microwell flat-bottom polystyrene plate, with RBD antigen	2-8 °C
C-2020	Dilution agent	1.8 g	1 bottle	Powder for preparing dilution solution	Room temperature
D-2020	Washing solution 10x	50 mL	1 bottle	Phosphate buffer at pH 7.4 with Tween 20 as surfactant	Room temperature
E-2020	Positive Control	12 µL	1 tube	SARS-CoV-2 anti-protein S antibody	2-8 °C
F-2020	Negative Control	12 µL	1 tube	Human IgG antibody not related to SARS-CoV-2	2-8 °C
G-2020	Secondary antibody	12 µL	1 tube	HRP-labeled anti-human IgG antibody	2-8 °C
H-2020	Substrate A for development	5 mL	1 bottle	Hydrogen peroxide solution	2-8 °C
I-2020	Substrate B for development	5 mL	1 bottle	3,3',5,5'-tetramethylbenzidine (TMB) in organic solvent	2-8 °C
J-2020	Stop solution	5 mL	1 bottle	Methane sulfonic acid solution	2-8 °C
N/A	Adhesive sealing for 96-well plates	1	1 film	Non-sterile adhesive plastic sealant	Room temperature

- 6 REQUIRED MATERIALS NOT INCLUDED IN THE KIT**
- Micropipette tips for 10, 200 and 1000 µL
 - Serological pipette for 5, 10 and 25 mL
 - 100 mL graduated cylinder
 - 1500 mL microcentrifuge tubes
 - 15- or 50-mL conic tubes
 - Reagent reservoirs
 - Deionized or distilled water
 - Microplate reader for reading at 450 nm

7 REAGENTS PREPARATION AND STORAGE
The kit is transported at room temperature (18-25°C) however, once the kit is in the laboratory it is suggested to refrigerate or separate the reagents according to their storage temperature. Prepare working solutions immediately before use. See reagent preparation section. Avoid using solutions with more than 24 hours of being prepared.

7.1 WASHING SOLUTION
Phosphate buffer at pH 7.4 with surfactant
Amount: 1x50 mL
Storage: Room temperature
Preparation: concentration 10X. Dilute content with 450 mL of deionized or distilled water, mixing properly before use.

7.2 DILUENT AGENT
3 % diluent agent in phosphate buffer with surfactant.
Amount: 1x1.8 g/60 mL
Storage: Room temperature no more than 25°C
Preparation: Transfer 60 mL of washing solution to the container marked as **DILUTION AGENT**, mix thoroughly before use. Once the diluting agent has been reconstituted it should be discarded at the end of the assay. It is not recommended to store it for use in further assays.

7.3 DILUENT SOLUTION OF ANTIBODIES
1% dilution agent in phosphate buffer with surfactant.
Amount: 1 x 12 mL
Storage: Room temperature no more than 25 °C
Preparation: Prepare a dilution 1:3 of **DILUTION AGENT** (prepared as in the step before) in a conical tube. Add 8 mL of washing solution in a conical tube and transfer 4 mL of dilution agent, mix thoroughly before use. It is not recommended to store it for use in further assays.

7.4 SECONDARY ANTIBODY
Anti-human IgG antibody coupled to HRP
Amount: 1x12 mL for 6 mL
Storage: 2-8°C
Preparation: Transfer the secondary antibody content to 6 mL of diluent solution of antibodies. Mixing before use. Once the antibody has been prepared it should be discarded at the end of the assay. It is not recommended to store it for subsequent use.

7.5 DEVELOPMENT SUSTRATES
Solution containing hydrogen peroxide (A) and 3,3',5,5'-tetramethylbenzidine (TMB) in organic solvent (B).
Amount: 2x5 mL
Storage: 2-8°C
Preparation: Mixing equal volumes of substrate A and substrate B for developing. The kit contains enough volume of substrate A and substrate B for a whole plate, if fewer wells are required, prepare the required volume only. **Prepare 15 minutes before use**, mix properly and allow to stand at room temperature. Once the developing substrate has been prepared, it should be protected from light and discarded at the end of the assay. It is not recommended to keep it for use in further assays.

7.6 STOP SOLUTION
Methane sulfonic acid solution
Amount: 1x5 mL
Storage: 2-8°C
Preparation: Ready-to-use

7.7 POSITIVE CONTROL
SARS-CoV-2 anti-protein S antibody
Amount: 12 mL for 0.3 mL
Storage: 2-8°C
Preparation: Transfer 0.3 mL of **DILUENT SOLUTION OF ANTIBODIES** to the tube labeled as positive control. Mixing thoroughly by gently pipetting before use. Once the antibody has been prepared it should be discarded at the end of the assay. It is not recommended to store it for use in further assays.

7.8 NEGATIVE CONTROL
Anti-human IgG antibody no related to SARS-CoV-2
Amount: 12 mL for 0.3 mL
Storage: 2-8°C
Preparation: Transfer 0.3 mL of **DILUENT SOLUTION OF ANTIBODIES** to the tube labeled as negative control. Mixing thoroughly by gently pipetting before use. Once the antibody has been prepared it should be discarded at the end of the assay. It is not recommended to store it for use in further assays.

8 SAFETY CONSIDERATIONS
Avoid skin to contact the reagents; wear gloves, laboratory coat and safety goggles when manipulating kit reagents. In case of accidental contact, wash thoroughly with water for 15 minutes at least. Due to potential risk of serum samples containing active SARS-CoV-2 viruses, heat treatment as described in sample preparation section is recommended.

9 SAMPLES
10 mL of human serum sample is required for determination in duplicate. Serums could be analyzed the same day as the sample collection, kept refrigerated (2 to 8°C) for up to 2 weeks (longer periods must be tested by the user) or keep it frozen (-20°C) for longer storage times. Analyzing hemolyzed samples is not recommended.

10 ESSAY PROCEDURE

10.1 Preparation of Regents
Prepare working solutions immediately before use. See reagents preparation section. Avoid using solutions with more than 24 hours of being prepared.

10.2 Preparation of Sample
a) Heating samples at 56°C in water bath is recommended. Such treatment is not necessary when samples are processed on a biosafety cabinet. No differences in assay performance have been observed between heat-treated and non-heat-treated samples.
b) Centrifugate to 2000 rpm for 3 minutes.
c) Dilute samples 1:1000, take 10 mL serum and add 900 mL of reconstituted **DILUTION AGENT**. Mixing properly before adding to the plate.

10.3 Essay Performance
a) Remove the 96-well plate from packaging and invert the excess of liquid on an absorbent material.
b) Add 100 mL of diluted samples 1:100 and the **CONTROLS** to each well, analyzing by duplicate all samples is recommended.
c) Incubate for 1.5 hours at room temperature.
d) Wash 6 times adding 200 mL of **WASHING SOLUTION** to each well and fully remove.
e) Add 50 mL of **SECONDARY ANTIBODY** previously diluted.
f) Incubate for 1 hour at room temperature.
g) Remove the **SECONDARY ANTIBODY**
h) Wash 6 times adding 200 mL of **WASHING SOLUTION** to each well and fully remove.
i) Bring to room temperature and prepare 15 minutes before use the **DEVELOPING SUBSTRATE**. Add 100 mL of **DEVELOPING SUBSTRATE** to each well.
j) Incubate at room temperature for 20 minutes avoiding light.
k) Add 50 mL/well of **STOP SOLUTION**.
l) Determine the OD via plate reader using a 450 nm filter and correction to 579 nm. Note: if the spectrometer does not have a filter to read at 570 nm, use another filter for correction (> 570 nm).

11 RECOMMENDATIONS FOR THE ASSAY
a) It is recommended to analyze all samples by duplicated, at least. The interpretation of results should be from the OD average of each sample.
b) Keep the **DEVELOPING SUBSTRATES** on the original container and avoid exposure to light.
c) A careful pipetting technique is necessary, and properly calibrated pipette devices are necessary to assure reproducibility of the test.
d) Different incubation times or temperatures to those stated in this manual could modify the results.
e) Avoid air bubble formation on wells, as this could result in a lower bind efficiency and a higher coefficient of variation. The reverse pipetting technique could be used for this purpose.
f) All reagents should be mixed thoroughly and gently before use. Avoid foaming.

12 QUALITY CONTROL
In order to assure validity of results, it is suggested each assay to include negative and positive internal controls, in addition to those contained in the kit. The OD average value of positive control should be not less than 1.5.

13 RESULTS INTERPRETATION
a) Calculate OD average value of controls and samples.
b) Determine the sample interpretation by comparing OD in the following table:

Interpretation	Interval	Results
Negative	Obtained value ≤ 0.5	The sample does not contain IgG antibodies related to SARS-CoV-2
Positive	Obtained value ≥ 0.6	The sample does contain IgG antibodies related to SARS-CoV-2
Limit	0.5 < Obtained value < 0.6	Repeat the test and interpret with other clinical tests

14 TEST LIMITATIONS
a) This test is for qualitative detection only. The test result should not be the only basis for clinical diagnosis and treatment. Confirmation of infection with SARS-CoV-2 should be combined with patient's clinical data and other diagnostic tests.
b) In the first week after the onset of SARS-CoV-2 infection, patient results could be negative for IgG, since it is an early stage for detection of these antibodies. In addition, patients with low immunity or other diseases affecting immune function, failure of main systemic organs, and usage of an inhibiting immune function drugs can also provide a negative result of IgG response to coronavirus. Previous infection with SARS-CoV or another coronavirus strain could provide a positive IgG response due to possible phylogenetic homologies among viral strains.
c) Bacterial or fungal contamination in serum samples or reagents, or cross contamination among reagents could provide incorrect results.
d) Deionized water treated with polyester resins could inactivate the HRP enzyme.

15 DEVELOPMENT FEATURES

15.1 Sensitivity
There are no international standard controls available to assess the immune response to SARS-CoV-2. To assess the

specificity of the assay, positive serums (from patients diagnosed as SARS-CoV-2 positive by qRT-PCR) and negative serums (obtained before pandemic and from SARS-CoV-2 negative patients by qRT-PCR) were analyzed to establish the OD endpoint and determine the presence of antibodies..

15.2 Accuracy
During the kit development, positive samples were evaluated in triplicate by two analysts independently. A coefficient of variation percentage (CV %) of OD less than 10% was obtained intra- and inter-assays.

15.3 Cross Reactivity
Cross reactivity with other diseases has not been determined.

16 REFERENCE
Cao X. COVID-19: immunopathology and its implications for therapy. Nat Rev Immunol. 2020;20:269–70.

Version 1.0; February 2021.